

HAWC Optical Calibration: 900 \rightarrow 300 Tanks

John A.J. Matthews and W. Miller (**UNM**)

Brenda Dingus (**LANL**)

Petra Huentemeyer (**MTU**)

johnm@phys.unm.edu

University of New Mexico

Albuquerque, NM 87131


Calibration system: *baseline design*



The HAWC array is now composed of 150 **tank-pairs**:


- As before we propose:
 1. two, identical laser sources: each source calibrates $\sim 1/2$ tanks
 2. lasers are coupled to optical fibers w/ beam expanders and 1:n fiber splitters
- What has changed is:
 1. one duplex optical fiber now goes to each **tank-pair** providing the calibration light to those tanks and 2-times the light path to monitor the end-to-end propagation time.
 2. at each **tank-pair** light is now split (only) 7-ways: 3 fibers to each tank and 1 fiber to provide the end-to-end time monitor.
 3. much shorter, $\sim 15\text{m}$ long, optical fibers go from the 1:7 splitters to a *diffuser* near each PMT.

Calibration system: *laser source (I)*



- What has changed is:
 1. time synchronization between the two (independent) laser sources now uses a PMT to define a t_0 ... that must be input to the HAWC DAQ.
 2. a beam shutter (in each laser source) ensures that calibration light comes from only one of the lasers (at any given time).
 3. a 1500 grit diffuser results in a more-uniform speckle pattern

Calibration system: *laser source (II)*



- Laser beam profile with 10x beam expander:
 1. a 1:91 optical fiber splitter (with $62.5\mu/125\mu$ optical fibers) occupies $\sim \pm 0.75\text{mm}$
 2. thus most (outgoing) signals should be within $\sim \pm 10\%$ in intensity.


Calibration system: *laser source (III)*

MacService - Mac Mini Repair - Mac Mini Hard Drives - Mac Mini Logic Boards 10/23/09 1:05 PM

MyService Repairs & Upgrades How It Works About Contact View Cart

Blog Business Services Education Program Support Your Account 877.MAC.DIRECT (622.3473)

Chat Now


Mac Mini
Repairs & Upgrades 

play video service details

Free Diagnostic
Included with all orders. Select this option if you're not sure which service you need. Free with all approved repairs. [read more](#)

Complete Diagnostic Free [Add to cart](#)
Includes free round trip shipping with all approved repairs.

Mac Mini Hard Drives
Store more music, photos, videos and files with a new hard drive upgrade. If your current hard drive is slow or clicking, a new hard drive upgrade can restore functionality to your mini. All MyService hard drives are new and carry full manufacturer's warranties.


 Includes installation of the new drive and transfer of your data. If your current drive is failing, a new operating system is installed. Data recovery options are available. Your old hard drive is returned to you. [read more](#)
[Play Video](#)

Intel Mac Mini Hard Drive Upgrades

160GB Hitachi 5400rpm	\$175	Add to cart
160GB Hitachi 7200rpm	\$205	Add to cart
250GB Fujitsu 5400rpm	\$225	Add to cart
250GB Hitachi 7200rpm	\$255	Add to cart
320GB Hitachi 7200rpm	\$295	Add to cart
500GB Samsung 5400rpm	\$325	Add to cart
500GB Seagate 7200rpm	\$345	Add to cart
120GB Solid State Drive	\$499	Add to cart
250GB Solid State Drive	\$899	Add to cart

- Possible calibration control computer:
 1. the Mac Mini can now be obtained with solid state disk drive
 2. a conventional Mac Mini computer (running UNIX) has been used to control one of the Auger FD relative calibration systems for ~ 18 months without problems.

Calibration system: *other studies*



- Mechanical adjustment of fiber intensity (or coupling into 1:n splitters):
 1. it is desirable for the return (end-to-end timing) signals to be relatively similar in amplitude.
 2. while fine control is difficult, mechanical separation of the fibers at a ST:ST coupler allows the transmitted intensity to be varied by $\sim 10\times$ per fiber ...

Calibration system: *continuing issues*

- What is the fiber to fiber uniformity and fiber (output) light intensity for **field distribution** options:
 - 1:7 optical fiber splitters, 62.5 μ m fibers, and one diffuser per PMT
 1. advantage is that the diffuser can be *as close as needed* to the PMT; the design goal is a maximum signal of $\sim 10^4$ PEs
 2. a possible dis-advantage is that there is one fiber per PMT VS one fiber per tank
 - 1:3 optical fiber splitters, 200 μ m fibers, and one diffuser per tank
 1. an advantage is that the larger fibers have smaller fiber to intensity variations
 2. the dis-advantage is that the diffuser must be $\gtrsim 2$ m from the PMTs; thus the PMT solid angle at the diffuser is small!
 3. it is unclear whether one diffuser/tank is a plus or minus in regards placement and tethering issues ...
- Or can we gain a significant ($\gg 2\times$) increase in intensity with a different laser?

Calibration system: *summary/conclusions*

- Continuing progress ...
- Nevertheless several ongoing studies ...
- Focus on few-tank **prototype array** is critical to help resolve the next set of issues ...